

Features	Dynamo SP10	Dynamo SP8	Dynamo SP6	Dynamo SP4
Non-linear fitting formulas	DSEsp, NAL, DSL	DSEsp, NAL, DSL	DSEsp, NAL, DSL	DSEsp, NAL, DSL
DSE personal profiles	5	3	3	1
Linear fitting formulas	DSElin, NAL-RP	DSElin, NAL-RP	DSElin, NAL-RP	DSElin, NAL-RP
Speech Rescue™	Yes	Yes	Yes	Yes
Speech Guard E	Yes	Yes	No	No
Single compression	No	No	Yes	Yes
Spatial Noise Management	Yes	No	No	No
Binaural Synchronisation (automatics)	Yes	Yes	No	No
Binaural Coordination (PB operations)	Yes	Yes	Yes	Yes
YouMatic	Premium	Advanced	Essential	No
Noise Management	Transient, Tri state	Tri state	Standard	Standard
Inium Sense feedback shieldSP	Yes	Yes	Yes	Yes
Free Focus	Premium	Advanced	Essential	Basic
Back dir	Yes	Yes	Yes	No
Power Bass	Yes	No	No	No
Music Widening	Yes	No	No	No
Voice Priorityi™	Yes	Yes	Yes	No
Special Purpose programs (music, lecture etc.)	Yes	Yes	Yes	No
Easy RECD™	Yes	Yes	Yes	No
Fitting Bands	9	8	6	4