

Irresistible Intuition

Call it a sixth sense. Or call it Enchant. The first Sonic instrument featuring the revolutionary SoundDNA platform, Enchant delivers natural sound in a brand new way.

We envision a world where all people can enjoy the sounds that enrich everyday life. With Sonic, everyday sounds better.



Zpower®

Change batteries just once per year with ZPower. Ask your Hearing Care Provider for details!

Bluetooth is a registered trademark owned by Bluetooth SIG, Inc. USA. Apple, the Apple logo, iPhone, iPad, iPod touch, and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

© 2018 Sonic Innovations, Inc. All rights reserved.

enchant™

Instinctively Smart



What makes Enchant so fascinating?

The first family of instruments built on the SoundDNA platform, Enchant brightens your everyday experiences with incredibly adaptive, flexible, and automatic technologies. Personalize Enchant with your own mix of features, accessories, and style.*

Capabilities That Captivate



SmartCompress

Sonic's renowned Speech Variable Processing (SVP) gets even smarter with SmartCompress technology. With enhanced processing power compared to our previous platform, SVP with SmartCompress accurately adjusts to your surroundings for stunningly natural listening experiences.



SPiN Management

The SPiN (Speech in Noise) Management system locates speech while suppressing background noise from various directions to help make conversations clear and comfortable.



Extended Dynamic Range

Extended Dynamic Range (EDR) helps keep speech and other sounds clear at movie theaters, performances, auditoriums, and other places where sudden dramatic sounds really make the moment.

Delight in Bits of Brilliance



Binaural Noise Management

This wireless option optimally reduces noise on either side of you, giving you a balanced sense of hearing in noisy environments.



Tinnitus SoundSupport

This feature helps provide relief from tinnitus by masking the ringing with a specialized signal, then increasing the volume of sounds you want to hear.



SmartMusic

Appreciate both the hushed harmonies and booming sounds of live performances. This program helps you enjoy undistorted music at its best.



Impressive IP Rating

Don't sweat the small stuff — like sweat itself. An IP68 rating means Enchant stands up to everyday exercise, spontaneous splashes, and general exposure to dust and debris.



Frequency Transfer

If you have difficulty hearing higher-frequency sounds, this feature transfers them to a lower region so you'll better hear subtle speech sounds like /s/ and /z/.



Telecoil

This handy feature picks up signals from your phone or from a loop system in venues such as theaters and places of worship that use this symbol.



Adaptive Feedback Canceller

This sophisticated system removes feedback signals — often before they're even heard.

*Not all features available on all technology levels or models. Ask your Hearing Care Provider for details.

Wireless That's Even More Wonderful




Dual-Radio System

The wow of wireless is even better in Enchant. Dual-Radio System uses 2.4 GHz wireless technology, so you don't need an intermediary device to stream sound from select digital sources to your hearing aids. And it uses **Near Field Magnetic Induction (NFMI)** technology, making communication between instruments faster.



Sonic Enchant and SoundLink 2 App are compatible with iPhone X, iPhone 8 Plus, iPhone 8, iPhone 7 Plus, iPhone 7, iPhone SE, iPhone 6s Plus, iPhone 6s, iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, 9.7-inch iPad Pro, 12.9-inch iPad Pro, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 4, iPad mini 3, iPad mini 2, iPad mini, and iPod touch (5th and 6th generation). Devices must be running iOS 9.3 or later. The app also supports Apple Watch. To download the Sonic SoundLink 2 App on iPad, go to the App Store, search for Sonic, and use the filter: iPhone only. Sonic SoundLink 2 App is compatible with devices powered by Android™ 6.0, Marshmallow or above. For information on compatibility, please visit www.sonici.global/compatibility.



SoundClip-A



RC-A



TV-A



Wireless Accessories

- » **SoundClip-A**
Enjoy stereo streaming with Bluetooth® compatible devices. It also functions as a remote microphone, remote control, and allows Enchant hearing aids to be used as wireless headphones for hands-free calls.
- » **SoundLink 2 App**
Change your listening program, adjust the volume, control tinnitus relief sounds, and switch between accessories, all from your iPhone® or Android™ smartphone.
- » **RC-A Remote Control**
Adjust volume and programs with a click. This remote is small enough for a pocket or purse.
- » **TV-A Adapter**
Stream favorite programs right through Enchant. This adapter easily connects to most TVs.

Enchant Style



Enchant miniRITE Enchant miniRITE T Enchant BTE 105




Invisible-in-canal Completely-in-canal In-the-canal In-the-ear

Enchant is available in five technology levels with premium features, plus simply brilliant models and colors.



beige taupe dark brown grey dark grey black* gold**

*Black color available in Enchant¹⁰⁰¹⁸⁰¹⁶⁰ only.
**Gold color available in the Enchant¹⁰⁰ miniRITE only.



black* beige light brown medium brown dark brown

*Black color available in IIC only.